



RENISHAW inVIA REFLEX MICRO-RAMAN

Instrument Quick Stats

- Accommodates solid, liquid and gaseous samples.
- Temperature controlled stage ranging from 25-1500°C.
- Computer-controlled switching between laser sources.
- Automatic alignment and optimization of laser source power.
- Laser spot size variable from 1-300µm.
- Facile switching between microscope viewing and Raman mode.
- Capable of confocal measurements with 2.5µm depth resolution.
- XYZ mapping sample stage with 0.1µm steps.
- High sensitivity, ultra-low noise CCD array detector.



Raman microscopy is a non-destructive analytical technique in which a material is excited by an incident laser at a known energy and subsequently characterized by plotting the intensity of scattered light versus the change in its energy. Every material has a characteristic Raman spectrum. A typical spectrum (for topaz) is shown above. Information about chemical composition, molecular structure and molecular interactions can be obtained using this technique.

The Renishaw inVia Micro-Raman is equipped with two laser excitation sources: a near infrared diode laser source (300mW) for excitation at 785nm and an Argon ion laser source (25mW) for excitation at 514nm. The instrument is designed for one-click switching between the sources, with computer-controlled reconfiguration and optimization. Calibration is automatic, using an internal reference sample. The micro-Raman is also equipped with a Leica microscope with binocular eyepieces and an integrated color video camera. The microscope is capable of making confocal measurements with 2.5µm depth resolution.

Contact Information

Thomas Wittberg
NEST Lab Manager • (937) 229-3001 • Thomas.Wittberg@udri.udayton.edu

Visit our website at <http://www.nestlaboratory.com>

NEST Laboratory • 300 College Park Drive • Dayton, OH 45469 • (937) 229 – 5776